



MyBreath iPhone App
Navigation and Flow Mockups

March 2, 2010
Version 2.0

Created by
Juan Juan




Outside of
iPhone App

Register


Main Menu

Appears after the splash page. Allows users to find the page they want to enter and an easy way to navigate to other pages. There will be a Menu button that will be in the top right corner on the other pages.



Info Page

Is linked to the I icon in the Main Menu. This page is for Company Info with a feedback link and website link at the end of page.




Breathe Page

A popup notice that will be displayed for a couple seconds.




Breathing Stage

No time limit in breathing. User can hear breathing, view moving breath graph and breaths per minute.



Special Effects
When user presses on Special Effects button a menu appears with items to choose from.



Menu
- Calm and Relax
- Ocean
- Cathedral

User chooses an item and menu disappears and special effects audio begins to play.


Artistic Rendering

Bar Graph is default and user can press button to choose different artistic display




Art Gallery

User can choose different artistic renderings from pre-installed list.




Kaleidoscope

User chooses Kaleidoscope and artwork is displayed in the iPhone display.




Breathing Stage

No time limit in breathing. User can hear breathing, view moving breath graph and breaths per minute.




Recording Stage

User presses Record Button and records Breathing for 1-5 minutes.



Save or Cancel


Cancel leads back to Breathe Stage and Save takes User to Saved Sessions and asked if they want to upload to website for full analysis.



Pace Stage


No time limit in breathing. User hears preinstalled pace, view moving breath graph and breaths per minute to match target breath per minutes.

User can control Pace by using arrows to increase or decrease breaths per minute.



Saved Sessions

After user saves breathing session. It is stored in this section. They can replay, delete, or upload to website. The date and time is displayed for when session was saved.



Register Page

After user uploads saved session to website. The Safari browser opens and user is asked to register or login to see full analysis of their breathing.

